# Avital 6276

**Blood Glucose Monitoring System**

**No Coding  
Required**

Model:6276

## USER GUIDE

11-62-0210-0077 V4-JAN22

## Explanation of symbols:

	Do not reuse
	Use by / Expiry date
	Lot number
	<i>In vitro</i> diagnostic medical device.
	Temperature limitation
	Caution (refer to accompanying documents). Please refer to safety-related notes in the manual accompanying this instrument.
	Consult instructions for use.
	Keep away from sunlight.
	Authorised representative in the european community.
	Manufacturer

### Important Notice to Reader of this Manual

This manual is written for the Avital Blood Glucose Monitoring System (6276).

**IVD** For self-use only.

Do not use this meter near cellular or cordless telephones in a call, walkie-talkies, garage door openers, radio transmitters, or other electrical or electronic equipment that are sources of electromagnetic radiation, as these may interfere with the proper operation of the meter.

# 1. Contents

<b>Section Topic</b>	<b>Page</b>
1. Table of Contents	1
2. Introduction	2-3
3. Appearance and Key Functions of the meter	4
4. Alternate Site Testing	5-7
5. Safety Information	8
6. Quick Testing Instructions	9

---

## **Preparing for Testing**

7. Installing the Battery	10
8. Settings (Date/Time/Alarm/ Memory Deletion /Voice volume)	11-14
9. Running a Control Solution Test	15-16

---

## **Blood Testing**

10. Running a Blood Glucose Test	17-19
11. Recalling the Memory and Viewing the Average	20-21
12. Expected Values Reference	22
13. Performance Evaluation Data	22
14. Transmitting Results	23

---

## **Maintenance**

15. Maintenance of Your System	24
16. Troubleshooting	24-25
17. Customer Service	26
18. Limitations	27-28

---

## **Important Additional Information**

19. Specifications	28-29
20. Ordering Supplies	29
21. Warranty Information	29

**NOTE:** For information about the name of the manufacturer of the lancing device and the lancets, please refer to the package inserts that came with your starter kits.

## **2. Introduction to the System**

### **Intended Use / Indications for Use**

The Avital Blood Glucose Monitoring System (6276) is intended for use in the quantitative measurement of glucose in fresh capillary whole blood from the finger and the palm, forearm, upper arm, calf or thigh. It is intended for use as an aid by either a healthcare professional or a patient diagnosed with diabetes mellitus for monitoring their blood glucose levels when self-testing. The Avital Blood Glucose Monitoring System is not intended for the diagnosis of or screening for diabetes mellitus, nor for use with neonates.

The alternative site testing (palm, forearm, upper arm, calf or thigh) in this system can only be used during steady-state blood glucose conditions.

### **Contents of the System**

This User Guide has been prepared to describe the proper use of the Avital Blood Glucose Monitoring System (6276).

Please read this User Guide and the package insert that comes with the Avital Glucose Test Strips before testing.

The system is available either as a meter alone or as a kit. If you have a meter only, you can purchase the system supplies from your system provider.



The Avital Blood Glucose Monitoring System (6276) contains the following items:

1. Avital 6276 Glucose Meter with 2 AAA batteries
2. User Guide
3. Carrying Case
4. Avital Glucose Test Strips
5. Lancets
6. Lancing Device
7. Package Inserts for Test Strips

You may purchase any control solutions that you need from your system provider.

The Avital Blood Glucose Meter uses Avital Blood Glucose Test Strips only. Neither the meter nor the test strips will work when used with any other brand of glucose products.

### 3. Appearance and Key functions of the meter

1. Test strip slot - When the strip is inserted into the slot, the meter will automatically turn on.
2. LCD Display - Guide you through the test using symbols and simple messages.
3. M Key - Power ON/OFF, also for memory recalling mode, please refer to manual for detailed function description.
4. C Key - Setting mode, please refer to manual for detailed function description.
5. Micro USB port - Cable connection and data transmitting.
6. Battery Compartment - Where batteries are located.
7. Ejector - Remove used strip.

\*If your meter has the back light feature, this will come on automatically during every operation to help user readability and to assist with the smooth operation of your meter.

#### The principle of the method

When glucose reacts with the reagents on the test strips, an electrical current is produced, which is proportional to the glucose concentration in the blood sample.

The glucose concentration is calculated by the meter and based on the current measured.

## 4. Alternate Site Testing (AST)

What is AST?

AST is the sampling from anatomical sites (parts of the body) other than the fingertip to check the blood glucose levels. This system allows you to test on the palm, forearm, upper arm, calf or thigh with the equivalent results to fingertip testing.

There are important limitations to AST. Please consult your healthcare professional before you use AST.

What is the advantage ?

Some patients may encounter discomfort when taking a blood sample from the fingertip as the finger tips have a high concentration of nerve endings. Other body sites do not have as many nerve endings, so you will not feel as much pain as at the fingertip.

When to use AST ?

Medication, stress, illness, food and exercise can affect blood glucose levels. Capillary whole blood at the fingertip can reflect test changes faster than capillary blood at other sites of body. If you test your blood glucose level during or immediately after a meal, physical exercise or stressful events, take the other sites of the body blood sample from your fingertip instead of from other sites.

Use AST only:


1. 2 hours or more after taking insulin;
2. 2 hours or more after a meal;
3. 2 hours or more after exercise.

Do not use AST if you are pregnant, or if you are aware that your glucose level is not as stable as usual, or if you think you have hypoglycemia (low blood sugar) or hyperglycemia (high blood sugar).


Do not use AST if you think your blood glucose is low and if your AST results do not match the way you feel.

## Testing Procedure for AST


1. Please unscrew the lancing device by turning the end cap counter clockwise. A clear cap will come with the package of a meter kit or a lancing device. The clear cap on the lancing device will make it easy for you to get a drop of blood for AST.


2. Insert a new lancet firmly into the lancet holder.

3. Twist off the protective top of the lancet.

4. Close the end cap of the lancing device and grasp the lancing device as shown. Slide into locking position. If necessary, please set the lancing device for a deeper puncture.

5. Select a soft, fleshy area on alternate testing sites that is clean and dry, away from bone and free of visible veins and hair.


6. Massage the selected area gently to increase blood flow to the puncture site. Clean the test location with an alcohol wipe or with soap and water.


7. Open the test strip vial. Take one test strip out of the vial and quickly cover the vial to keep air out.

8. Insert a test strip into the meter while the meter is off. The meter will turn on automatically. A blinking blood icon will be displayed on the screen. If no action is taken in 2 minutes, the meter will turn off automatically.

9. Press and hold the loaded lancing device against the other test sites for a few seconds, then press and release the button.

## 5. Safety Information

Please use this device only for the intended use described in this user guide.

Before using this system to test your blood glucose please read the instructions carefully and practice the test when you first use this system. Do a quality check on the system following the instructions and consult with your healthcare professionals for questions or problems.

Be aware of the safety of young children or handicapped persons near you when you conduct a glucose test using this system.

Please keep the test strip vial away from children.  
The test strips and vial cap can present a choking hazard.


Please be cautious when removing the lancet. Take the lancet out carefully. Always place the protective cap back on the exposed tip or (if using a twist lancet) place the cap on a hard surface and push the exposed tip into the protective cap and discard of properly.

Never try to disassemble the meter in any circumstance. If your meter is not working properly, or if you need the Technical Support, please contact your local distributor for help.


## 6. Quick Testing Instructions

(For detailed instructions start on page 15 for Control Solutions and page 17 for Blood Testing.)


1. With the meter turned off, insert a new test strip into the meter. A blinking blood icon will be displayed on the screen.


2. Lance the finger and let a blood drop form.


3. Apply the blood drop to the front edge of the test strip when the blood drop icon is still blinking on the meter. Wait for 6 seconds, and the meter will display the test result.


4. Remove the used strips by hand or by pushing the ejector and the meter will power off.


## 7. Installing the battery


1. Turn off the meter by pressing the "M" button and hold for 2 seconds before replacing the batteries. Slide down the battery cover on the back of the meter by pushing the cover in the direction of the arrow and put the cover aside.

2. Remove the two used batteries.

3. Insert 2 new batteries with correct +/- direction. (Battery: AAA 1.5V 2 Alkaline 24A LR03)


4. Put the battery cover back in place by pushing the cover in the opposite direction of the arrow to close the cover into position.

### Note:

1. It is suggested that batteries need to be replaced when the battery icon appears empty. If you keep using the meter, the battery icon will start flashing and an error message E\_3 will appear on the screen.
2. Note that replacing the battery will not reset the stored testing results.
3. You need to reset the time and date after the batteries are replaced.
4. If there is any unexpected symbol on the display, please follow the above procedure to replace the batteries again.
5. Please dispose of batteries according to your local ordinances.

## 8. Settings (Date / Time / Alarm / Memory Deletion / Voice Volume)

The LCD screen will show as below:


1. Appears with the test result in mg/dL or in mmol/L.
2. Appears when you recall the memory.
3. Appears when you are in setting mode.
4. Appears for control solution test flag (cotr).
5. Indicates the meter is ready to take the blood sample when it flashes.
6. Indicates the battery status.
7. Indicates current displayed result is an average.
8. Appears when alarm is on.
9. Month
10. Day
11. Hour
12. Minute


1. Set the Year - To go into the "Setting Mode" you must press and hold down the M Button. Then release the M Button and the first number on the screen will be the year. Press C key to adjust the year until the desired year is displayed and then press the M key to confirm the year setting. When pressing the C key, the year will be adjusted in the range of 2014 to 2029. After the year setting is confirmed, you will see the month segment flashing on the screen.


2. Set the Month - Press C key to adjust the month until the desired month is displayed and then press the M key to confirm the month setting. You will then see the day segment is flashing on the screen.


3. Set the Day - Press C key to adjust the day until the desired day appears and then press the M key to confirm the day setting. You will then see the hour segment is flashing on the LCD screen.


4. Set the Hour - Press C key to adjust the hour until the desired hour appears on the screen and then press the M key to confirm the setting. You will then see the minute segment is flashing on the screen.

5. Set the Minute - Press C key to adjust the minute until the desired minute appears on the screen and then press the M key to confirm the setting. You will then see an "AL" displayed on the screen center.


6. Set the Alarm Reminder - Press C key to toggle between the alarm reminder setting or no-setting mode. If you see the LCD display a small flashing icon of alarm clock in the lower left corner with a time segment display, it means that you are in the alarm setting mode.

If you see only AL without the icon of the small alarm clock, it means that you are not in the alarm setting mode. Press M key to confirm the alarm setting mode. If you select the alarm setting mode, alarm clock icon will stop flashing and the hour of the time segment will start flashing. Press C key until the desired alarm hour appears and then press M key to confirm the hour setting. Then the minute segment will start flashing. Press C key until the desired minute appears on the screen and then press M key to confirm the minute setting. You will see the measuring unit flashing on the upper left corner of the LCD display. Each time you press the C key to adjust the minute.

7. Delete Memory - Press C key to toggle between a flashing "dEL" for memory deletion or a non-flashing "dEL". If you want to delete all memory, press M key while the "dEL" is flashing. If you do not want to delete all memory, press M key when the "dEL" is not flashing.

7

**Note:**


1. When you change the date and time backward, you will not change the test results in the meter memory.
2. You need to move through the language, year, month, day, hour, minutes, alarm, measurement unit, dEL, to save the meter settings and turn off the meter.
3. The average readings in the meter memory are calculated from the results obtained during the 7, 14, 28, 60, and 90 calendar days preceding the current date and time settings.
4. Your meter can show results either in "mg/dL" or "mmol/L" units. It is very important to use the correct unit of measurement to properly manage your diabetes. If you live in the Europe, you should use mmol/L; your results will always have a decimal point. The unit of measurement will be set to mmol/L by default. If you do not see a decimal point in your test results, it indicates the measuring unit is in mg/dL.

## 9. Running a Control Solution Test


### IMPORTANT:

Always make sure you press C key for Control Solution test, otherwise the control solution test result will be stored in the memory, and affect your “AVG” results.


1. Insert a strip into the meter. A blinking blood icon will be displayed on the screen.

2. Press the C key, and the meter will display a letter “cotr” on the down right corner of screen, indicating the control test mode. The test result will not be stored in the memory.

3. Squeeze a small amount of Control Solution on a flat surface and apply to the front edge of the test strip.

4. Wait for 6 seconds, and the meter will display the result.

5. Remove the used strips by hand or by pushing the ejector and the meter will power off.


 **IMPORTANT!!!**


1. You must press the C key to distinguish the control solution test from the blood glucose test. Do not perform the blood test in Control Solution mode. (If you perform the blood test in Control Solution mode, the test result will not be stored in the memory.)
2. Use only Avital Control Solution with the Avital Blood Glucose Meter.
3. Avital Control Solutions are used to check that the meter and the test strips are working together as a system and that you are performing the test correctly.
4. Avital Control Solution contains a known amount of glucose that reacts with test strips. The Normal level Control Solution are intended to check the monitoring system in different measurement ranges.
5. Shake the vial, discard the first drop of control solution, and wipe off the dispenser tip to ensure a good sample and an accurate result.
6. Use only for 6 months after first opening. Record the open date on the Control Solution vial. Discard after 6 months.
7. Compare your Control Solution test results with the expected range printed on the test strip vial label. If your glucose control results fall outside the expected range; repeat the test. Results that repeatedly fall outside the expected range may indicate:
  - (1) You may not be doing the test correctly; repeat the test, by shaking the control solution vial well and carefully following instructions. Please make sure the test is done within the temperature range.
  - (2) The Control Solution is expired or contaminated.
  - (3) The test strips are damaged or outdated.
  - (4) You may have applied the Control Solution before the blood drop appears on the screen. This will cause incorrect glucose measurement.
  - (5) Meter malfunction.
8. The recommend temperature range for control solution test is 50-86°F (10-30°C). Always allow at least thirty (30) minutes to let meter, test strip and control solution adjust to the stable temperature before performing a control test.

## 10. Running a Blood Glucose Test


1. Please unscrew the lancing device by turning the end cap counter clockwise.


2. Insert a new lancet firmly into the lancet holder.

3. Twist off the protective cover of the lancet.

4. Replace and twist the adjustable depth cap onto the body of the lancing device until you hear the click. When the cap is properly closed in locked position, you should see the arrow and release button are aligned. (see Figure 3-4). Holding the cap and select lancing depth.

Slide the lancing device into locking position: gently pull the gray sliding barrel away from the top of the lancing device until a soft click is heard.


5. Wash your hands in warm, soapy water. Be sure to rinse and dry well.


6. Open the test strip vial. Take one new test strip out of the vial and recap the vial quickly and firmly.


7. While the meter is off, insert a new test strip into the meter's strip slot with the contact bars facing upwards. The meter will turn on with a full screen display. A blinking blood icon will be displayed on the screen. If no action is taken in 2 minutes, the meter will turn off automatically.


8. Use the adjustable lancing device or the disposable lancet to take a blood sample from the fingertip and let it form a small round blood drop. To avoid body fluid contamination, it is recommended that you do not use the first blood drop for glucose testing.

9. Apply a small blood drop to the front edge of the test strip, and the blood should be pulled into the confirmation window before the meter begins to count down. Do not push your finger against the test strip or try to apply a smeared blood sample.

10. The meter will start counting down for 6 seconds and your test result will be displayed on the screen. Remove the used test strip by hand or by pushing the ejector button. The meter will turn off and the test result is stored automatically.

**Important:**

1. Never reuse a test strip or a lancet. An "E\_5" error message will be displayed in the center of screen if a used strip is inserted into the meter.
2. Discard the test strip and lancet immediately in a puncture-proof container with a lid after use.

## **11. Recalling the Memory and Viewing the Average**

### **A. Using the Meter Memory**

Your blood sample test result will be automatically stored in the meter's memory. The meter can store up to 500 of the most recent test results. You can also view the average of test results for various periods of 7, 14, 28, 60 and 90 days. The control solution test results will not be stored in the memory. (Note: Please make sure that you press C key when you run a Control Solution test, so that the control test result will not be added to the memory.)

### **B. Recalling the Test Results from the Meter Memory**

1. With meter off, press M key to turn on the meter. After a brief full screen display, the meter will display the total number of test data stored in memory. Then the screen will display "001" as the sequential number of the most recent test result in memory. The memorized test data will be displayed next. If there is no memorized result in the meter, "---" is displayed. The meter will turn off after 60 seconds if no action is taken.
2. The screen will display each of the memorized test data in sequence. Press M and C key to review your test results in the memory forwards and backwards.
3. To turn off the meter, press the M key for 2 seconds, or the meter will turn off automatically after 60 seconds.

### C. Viewing the Average (AVG) of Test Data in Memory

1. To display the average (AVG) test result, please turn the meter off. Press C and M keys at the same time and hold for about 2 seconds, wait for the AVG to be displayed on the screen, then release the two buttons. The meter will be in the average mode with “AVG” font displayed.
2. When the “AVG” is flashing at the lower left corner of the screen, the number 7 will be displayed under AVG and the average of the last 7 days test results will be displayed in the center of the screen. The number of tests completed in the last 7 days will also be displayed in the lower right hand corner of the screen. If you want to see the average of 14, 28, 60 and 90 days, keep pressing the M key to move the LCD display forward from 7 day average to 90 day average. The average display will move from 90 days backward to 7 days if you continue to press C key.
3. Press the M key and hold for only 2 seconds to exit the recall mode and turn off the meter, or the meter will automatically power off after 60 seconds.
4. When “---” is displayed, showing that there are no test results in the memory.

### Calculations:

The 7 days average is the average of the last 7 day test results.

The 14 days average is the average of last 14 day test results.

The 28 days average is the average of last 28 day test results.

The 60 days average is the average of last 60 day test results.

The 90 days average is the average of last 90 day test results.

### Important:

1. Never reuse a test strip or a lancet. An “E\_5” error message will be displayed in the center of screen if a used strip is inserted into the meter.
2. Discard the test strip and lancet immediately in a puncture-proof container with a lid after use.

## 12. Expected Values Reference<sup>1</sup>

Blood glucose levels normally will vary from time to time depending on food intake, medication dosages, health, stress or exercise. Consult your physician or healthcare professional for the target glucose value appropriate for you.

Expected plasma blood glucose values for normal, nondiabetic adults are as follows.

Before eating < 5.6 mmol/L

Two hours after meals < 7.8 mmol/L

Avital 6276 meter gives plasma equivalent results.

## 13. Performance Evaluation Data

### Precision

Standard deviation (SD) for each glucose concentration < 100 mg/dL (5.55 mmol/L) and coefficient of variation (CV) for each glucose concentration  $\geq$  100 mg/dL (5.55 mmol/L) is  $\leq$  5.0 mg/dL (0.278 mmol/L) and  $\leq$  5.0%, respectively.

Intermediate precision

Control Solution Level (mg/dL)	Low (30-50)	Normal (96-144)	High (280-420)
Pooled Mean (mg/dL)	44.3	100.3	350
SD	3.0	3.2	8.5
CV (%)	—	3.2%	2.4%

### Repeatability

Blood Glucose (mg/dL)	30-50	51-110	111-150	151-250	251-400
Pooled Mean (mg/dL)	44.4	80.6	132	198.9	309.7
SD	3.0	2.3	4.1	5.8	9.0
CV (%)	—	2.9%	3.1%	2.9%	2.9%

### System Accuracy

For glucose concentration <100 mg/dL (5.55 mmol/L)		
Within $\pm$ 5 mg/dL (Within $\pm$ 0.28mmol/L)	Within $\pm$ 10 mg/dL (Within $\pm$ 0.56mmol/L)	Within $\pm$ 15 mg/dL (Within $\pm$ 0.83mmol/L)
119/180(66.1%)	163/180(90.6%)	180/180(100%)

For glucose concentration $\geq$ 100 mg/dL (5.55 mmol/L)		
Within $\pm$ 5%	Within $\pm$ 10%	Within $\pm$ 15%
233/420(55.5%)	387/420(92.1%)	417/420(99.3%)

For glucose concentrations between 44.7 mg/dL (2.48 mmol/L) and 459 mg/dL (25.5 mmol/L)
Within $\pm 15$ mg/dL (0.83 mmol/L) or $\pm 15\%$
597/600 (99.5%)

The Avital Blood Glucose Monitoring System meets the requirements for System Accuracy as stated in EN ISO 15197:2015.

### User Performance

A study evaluating glucose values from fingertip, palm, forearm, upper arm, calf, and thigh capillary blood samples obtained by 100 lay persons showed the following results: Fingertip 100% / palm 100% / forearm 100% / upperarm 100% / calf 100% / thigh 100% within  $\pm 15$  mg/dL ( $\pm 0.83$  mmol/L) of the medical laboratory values at glucose concentrations below 100 mg/dL (5.55 mmol/L), and fingertip 100% / palm 98.6% / forearm 100% / upperarm 100% / calf 98.6% / thigh 100% within  $\pm 15\%$  of the medical laboratory values at glucose concentrations at or above 100mg/dL (5.55 mmol/L).

## 14. Transmitting Results

### Note:

Avital 6276 Glucose Meter allows you to transfer the test results stored in its memory to your personal computer. However, you will need to order software and a Data Download Cable separately from your distributor. The software and its data download procedure available at distributor and on internet. You can also check with your distributor for the software download option. The meter still keeps the results in the memory after transmitting.

1. With the meter off, hold the M and C keys down at the same time for 5 seconds until "PC" is displayed.
2. Press the C key to start transmitting. The meter will display a flashing "PC" during the transmission process. When the "PC" on the screen stops flashing, it indicates the transmission is completed. The meter will be power off if no action is taken in 60 seconds or M key is pressed for 2 seconds.

## 15. Maintenance of your system

Please, treat this meter with proper care, and keep it in good condition.

1. Store your meter in the protective case, in a clean dry place from 50°F (10°C) to 104°F (40°C).
2. Always clean your meter after use. Wipe and clean the surface of the meter with a soft cloth that has been slightly dampened with mild detergent.
3. Please handle with care and do not drop the meter.
4. If this meter is used by healthcare professional infection control policies should be strictly followed.

We suggest you should periodically compare the test system to another test system which is well maintained and monitored by a healthcare provider

## 16. Troubleshooting :

The following is a summary of all the Error Messages.

These messages help to identify certain problems, but do not appear in all cases when a problem has occurred.

Improper use may cause an inaccurate result without producing an error message or a symbol. In the event of a problem, refer to the information in the table under Solution.

### Error messages:

**E\_1:** The temperature is too low.

**E\_2:** The temperature is too high.

**E\_3:** Battery Low.

**E\_4:** Memory damaged.

**E\_5:** The strip is wet or used.

**E\_6:** Error in meter or strip.

**E\_7:** The blood sample is not sufficient.

**HI :** The glucose level is too high.

**LO:** The glucose level is too low.

<b>Problem</b>	<b>Possible Cause</b>	<b>Solution</b>
(1) The meter displays an E_1 error message.	The meter is operating in an ambient temperature below 10 °C or 50°F which is the lowest range of operating temperature.	Repeat the test after the meter and strip are placed in a warmer environment and allow the meter to warm up for a while before retesting.
(2) The meter displays an E_2 error message.	The meter is operating in an ambient temperature above 40 °C or 104°F which is the cap of operating temperature.	Repeat the test after the meter and strip are placed in a cooler environment and allow the meter to cool down for a period before retesting.
(3) The meter displays an E_3 error message.	The battery is too low to operate the meter.	Please replace the two AAA batteries.
(4) The meter displays an E_4 error message.	The memory chip of the meter could be damaged or malfunctioning.	The meter can perform the glucose measuring without storing the test results into the meter's memory. Please write down the test results into the data logger before you call a customer service representative.
(5) The meter displays an E_5 error message.	1.If you apply the sample before blood drop icon appears on the screen. 2.The inserted test strip has been wet or used.	Please check the test strip to see if it is damaged or used. In either case, please discard the strip and repeat the test using a new strip.
(6) The meter displays an E_6 error message.	The meter is not working properly, either because of a defective meter or a defective strip.	Repeat the test with a new test strip. If E_6 continues to show up on screen, please call a customer service representative.
(7) The meter displays an E_7 error message.	The blood sample is not sufficient.	Please re-test by inserting a new test strip ensuring an adequate amount of blood comes into contact with the test strip.
(8) The meter displays HI on screen.	The test result is higher than the measuring range of 33.3 mmol/L.	Test again following the user guide for correct glucose measurement process. If you see HI again, please call your medical doctor for advice immediately.
(9) The meter displays LO on the screen.	The test result is lower than the measuring range of 1.1 mmol/L.	Test again following the user guide for correct glucose measurement process. If you see LO again, please call your medical doctor for advice immediately.

## 17. Customer Service

If you need assistance with your Avital Blood Glucose Monitoring System, please contact:

**ETTM LTD.**

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For questions related to your health condition,  
please call your doctor or healthcare professional

 **OK Biotech Co., Ltd.**

No. 91, Sec. 2, Gongdao 5th Road, 30070 Hsinchu City, TAIWAN

 **MDSS GmbH**

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30175 Hannover, GERMANY

  
0123

## 18. Limitations

1. Do not use for neonatal blood glucose testing.
2. Hematocrit in the range of 20% to 60 % does not affect the glucose results. If you do not know your hematocrit level, consult your healthcare professional.
3. Hemoglobin levels of 3.13 mmol/L (20000 mg/dL) did not interfere with blood glucose test results.
4. Cholesterol levels up to 12.9 mmol/L (500 mg/dL) and Triglycerides up to 99.3 mmol/L(3000 mg/dL) has been shown not to affect glucose results.
5. Interference was observed for therapeutic levels of Glutathione and Pralidoxime Iodide. No interference was shown for Uric acid, L-Dopa, Acetaminophen, Ascorbic acid and Gentisic acid in normal therapeutic levels. However, higher concentrations in blood may cause incorrect results:

Chemical	Limitation	Chemical	Limitation
1. Acetaminophen	≤ 8.0 mg/dL	14. Hydroxyurea	≤3.0 mg/dL
2. Ascorbic acid	≤ 5.0 mg/dL	15. Ibuprofen	≤50 mg/dL
3. Aspirin	≤ 120 mg/dL	16. Icodextrin	≤13 mg/dL
4. Bilirubin	≤ 90 mg/dL	17. L-dopa	≤ 10 mg/dL
5. Cholesterol	≤ 500 mg/dL	18. Maltose	≤ 900 mg/dL
6. Creatinine	≤ 6.3 mg/dL	19. Methyl dopa	≤ 3.0 mg/dL
7. Dopamine	≤ 2.0 mg/dL	20. Pralidoxime Iodide	≤25 mg/dL
8. EDTA	≤800 mg/dL	21. Salicylate	≤120 mg/dL
9. Galactose	≤ 900 mg/dL	22. Tolazamide	≤100 mg/dL
10. Gentisic acid	≤ 5.0 mg/dL	23. Tolbutamide	≤ 400 mg/dL
11. Glutathione	≤53 mg/dL	24. Triglyceride	≤ 3,000 mg/dL
12. Haemoglobin	≤20,000 mg/dL	25. Uric acid	≤ 11 mg/dL
13. Heparin	≤8000 U/dL	26. Xylose	≤100 mg/dL

6. The test strips may be used at altitudes up to 11,161 feet (3,402m) without an effect on test results.
7. Persons suffering from severe dehydration should not be tested using a capillary whole blood sample.
8. Test results below 3.3 mmol/L(60 mg/dL) indicate low blood glucose (hypoglycemia). Test results greater than 13.3 mmol/L (240 mg/dL) indicate high blood glucose

(hyperglycemia). If you get results below 3.3 mmol/L (60 mg/dL) or above 13.3 mmol/L (240 mg/dL), repeat the test, and if the results are still below 3.3 mmol/L (60 mg/dL) or above 13.3 mmol/L (240 mg/dL), please consult your healthcare professional immediately<sup>2</sup>.

9. Inaccurate results may occur for individuals experiencing a hyperglycemic-hyperosmolar state, with or without ketosis. Critically ill patients should not be tested with a blood glucose meter.

## 19. Specifications

Test: Glucose

Sample: Whole Blood

Test sites: Fingertip, palm, forearm, upper arm, calf or thigh

Measurement time: 6 seconds

Code required: No code device

Measurement Range: 1.1 to 33.3mmol/L.

Batteries: 2 x AAA 1.5 V Alkaline 24A LR03

Operating Temperature: 10°C to 40°C (50°F to 104°F)

Humidity: 10% to 85% RH

Width: 48 mm

Length: 99 mm

Thickness: 14 mm

Weight: ~72g

Memory: 500 measurement results with date and time

Auto power off after 2 minutes without action when a strip is inserted, or 60 seconds if there is no strip inserted.

Traceability of glucose monitoring system: The accuracy of Avital Blood Glucose Monitoring System was assessed by comparing blood results obtained by patients with those obtained using a YSI Model 2300 Glucose Analyzer, which is calibrated by YSI calibrator solution, refer to NIST SRM (standard reference material) # 917A Clinical Dextrose.

The device has certified to meet the following standard:

IVDD 98/79/EC, EN 61010-1, EN 61010-2-101, EN 61326-1, EN 61326-2-6, ISO 15197:2013

(P.S. Please refer to the package insert that came with your test strips, and you will also find information on how the system works and on the performance characteristic reference.)

## **20. Ordering Supplies**

To order supplies, please contact the Distributor:

**ETTM LTD.**

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## **21. Warranty Information**

Your Avital 6276 Glucose Meter is warranted to be free of defects in materials and workmanship for one year from the date of the original purchase. This warranty does not cover device failure due to owner misuse or negligence, or normal wear and tear.

If you have a question about your Avital 6276 Glucose Meter or this warranty, please contact:

**ETTM LTD.**

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Lancets and lancing device are medical devices under the Directive 93/42/EEC.

Manufacturer (Lancets) :

Please refer to lancet packaging for the details.

### **Reference:**

1. American Diabetes Care, January 2015, vol. 38 no. Supplement 1 S8-S16.

2. American Diabetes Association website:

<http://www.diabetes.org/living-with-diabetes/treatment-and-care/blood-glucose-control/hyperglycemia.html>

<http://www.diabetes.org/living-with-diabetes/treatment-and-care/blood-glucose-control/hypoglycemia-low-blood.html>

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